

SIGNAL COUNTING FOR IN SITU HYBRIDIZATION

ABSTRACT

Fluorescently tagged nucleic acid probe signals are counted in
5 biological specimens by determining a ratio of signals from a test probe to
signals of a reference probe. Probe signals need not be counted with
reference to cells, nuclei, or nuclear contours. Gene amplification or deletion
can thus be detected by analyzing the ratio. Successive image slices are
obtained by confocal microscopy, and the images are digitized. The digital
10 images are transformed and analyzed to combine contiguous fluorescent
signal segments in successive optical sections to identify discrete probe
signals, or spots. Spots overlapping in the axial and transverse dimensions of
a three-dimensional representation of the biological specimens can be
distinguished. A graphical user interface presents various features for
15 consideration by a user, who can provide guidance to a computer system
counting the spots. Various features directed to identifying spot clusters and
autofluorescent material can increase accuracy of spot counting.